

MaxPlus Transfusion Cooler

Designed and validated for
hospital transfusion services &
bedside storage of Red Blood
Cells or Thawed Plasma

CAPACITY

1 to 10 units



TEMPERATURE

1-6°C



QUALIFICATION

Up to 24 hrs



Reliable Blood Storage & Transport

Spend less time packing. Move blood faster.

Purpose-built pack-outs eliminate bench-time conditioning and special handling. With minimal training required, teams can pack quickly and consistently improving response times while reducing operational complexity. Each cooler comes pre-configured with the appropriate coolant packs, so it works reliably every time.


Built for Daily Blood Bank Use

Durable, reusable, and made to last.

Constructed from rugged Expanded Polypropylene (EPP), the MaxPlus Transfusion Cooler is engineered for long service life under heavy clinical use. It is chemically durable, and easy to clean using standard lab and OR wipes without degrading performance or structural integrity.

At-a-Glance

- Proven in field
- Reliable temperature protection
- Built to last
- Rigorously tested and qualified
- Pack and issue blood units in 30 seconds
- Designed for blood bankers



“As we moved the outpatient transfusion center across the lot, we realized that we needed a new way to transport blood products. This is when we ordered MaxQ Coolers. They were a better product than I ever knew existed. We are also using MaxQ’s frozen specimen pouches to pick up frozen samples from our numerous off-site locations.”

**Blood Bank Supervisor
Stillwater Medical Center**

MaxPlus Transfusion Cooler

Ergonomic by Design

Move blood safely, comfortably, and without distraction.

Designed for real hospital workflows, the MaxPlus Transfusion Cooler features smooth, rounded edges with no sharp corners, molded-in ergonomic handles, and optional shoulder straps for hands-free transport. The system also supports wheel attachment for longer distances and high-volume movement. Every detail is engineered so staff can focus on delivering the right blood product at the right time, helping save lives when seconds matter.



Product	Red Blood Cells / Liquid Plasma Units
Capacity	1-10 Units
Temperature	1-6°C
Pre-Qualification	Up to 24 hours
Coolants Needed	MTC13-GK (5 x coolant bottles)

CONTACT US



sales@flymaxq.com



405-466-5629



packmaxq.com/contact-us